

## Recombinant Human HSP60

|                           |   |
|---------------------------|---|
| <b>Catalogue No.:</b>     | P0250   |
| <b>Synonyms:</b>          | 60 kDa chaperonin, Chaperonin 60, cpn60, GROEL, heat shock 60kD protein 1 (chaperonin), heat shock 60kDa protein 1 (chaperonin), Heat shock protein 60, heat shock protein 65, HLD4, HSP60, HSP-60, HSP60SPG13, HSP65, HSPD1, HuCHA60, Mitochondrial matrix protein P1,60 kDa heat shock protein, mitochondrial, P60 lymphocyte protein, short heat shock protein 60 Hsp60s1, spastic paraplegia 13 (autosomal dominant), SPG13 |
| <b>Species:</b>           | Human   |
| <b>Uniprot ID:</b>        | P10809  |
| <b>Expression Region:</b> | 225-573   |
| <b>Source:</b>            | E.Coli  |
| <b>Tags:</b>              | GST tag   |
| <b>Molecular Weight:</b>  | 65.3 KDa  |
| <b>Purity:</b>            | Greater than 95% by SDS-PAGE gel analyses   |
| <b>Formulation:</b>       | Lyophilized from a 0.2um filtered solution in PBS with 5% trehalose, pH7.4  |
| <b>Reconstitution:</b>    | Sterile distilled water   |
| <b>Storage:</b>           | -20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile conditions after reconstitution   |


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)

**Image:****Safety note:**

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)