

## **Recombinant Human PRDX4**

Catalogue No.: P0798

**Synonyms:** Antioxidant enzyme 372, Antioxidant enzyme AOE372, AOE37 2,

AOE37-2, AOE372, EC 1.11.1.15, Peroxiredoxin IV, Peroxiredoxin-4,

Peroxiredoxin4, PRDX 4, Prdx4, PRDX4\_HUMAN, PRX 4, Prx IV, Prx-IV,

PRX4, PrxIV, Thioredoxin dependent peroxide reductase A0372,

Thioredoxin peroxidase

**Species:** Human

Uniprot ID: Q13162

**Expression Region:** 38-123

**Host:** E.Coli

**Tags:** N-terminal His-IF2DI Tag

Molecular Weight: 30.7 kDa

**Purity:** Greater than 95% by SDS-PAGE gel analyses

Formulation: Lyophilized from a 0.2 μm filtered solution in 10 mM Hepes, 500 mM NaCl

with 5% trehalose and 0.06% proclin, pH7.4

**Reconstitution:** Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 μg/ml

in sterile distilled water

**Storage:** -20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile

conditions after reconstitution

## Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

## **FineTest Proteins**



Image:



Safety note:

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

Tel: (0086)027-87384275 Fax: (0086)027-87800889 **www.fn-test.com**