

Recombinant Human CGA


Catalogue No.:	P1313
Synonyms:	Anterior pituitary glycoprotein hormones common subunit alpha, CG
	alpha, CG-alpha, cga, Choriogonadotropin alpha chain, Chorionic
	Gonadotrophin Alpha, Chorionic Gonadotrophin Alpha Subunit, Chorionic
	gonadotrophin subunit alpha, chorionic gonadotropin
Species:	Human
Uniprot ID:	P01215
Expression Region:	25-116
Host:	E.Coli
Tags:	N-terminal His-IF2DI Tag
Tags: Molecular Weight:	N-terminal His-IF2DI Tag 30.8 KDa
-	
Molecular Weight:	30.8 KDa
Molecular Weight: Purity:	30.8 KDa Greater than 95% by SDS-PAGE gel analyses
Molecular Weight: Purity:	30.8 KDa Greater than 95% by SDS-PAGE gel analyses Lyophilized from a 0.2 μm filtered solution in 10 mM Hepes, 500 mM NaCl
Molecular Weight: Purity: Formulation:	30.8 KDa Greater than 95% by SDS-PAGE gel analyses Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose and 0.06% proclin, pH7.4
Molecular Weight: Purity: Formulation:	30.8 KDa Greater than 95% by SDS-PAGE gel analyses Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose and 0.06% proclin, pH7.4 Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei,
China (430206)Tel: (0086)027-87384275Fax: (0086)027-87800889www.fn-test.com



Image:



Safety note:This product is intended for research and manufacturing uses only. It is not
a diagnostic device. Product degradation will result from multiple
freeze/thaw cycles. It is suggested that the antigen be stored in use size
aliquots and thawed just prior to use. This material has been inactivated,
however as with all biological materials, it should be handled as potentially
infectious. The user assumes all responsibility for care, custody and control
of the material, including its disposal, in accordance with all regulations.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei,
China (430206)Tel: (0086)027-87384275Fax: (0086)027-87800889www.fn-test.com